




RDE measuring head for Professional VA/CVS instruments with high chemical resistance

6.1256.040

Measuring head for operation with rotating disk electrodes. Solvent-resistant version for applications in nonaqueous media.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.1256.040


Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.1244.050

Drive belts for rotating disk electrodes

For Professional VA/CVS instruments (3 pcs.).



1 PCS

6.1829.030


PTFE tubing / M6 / 21 cm, tapered



1 PCS

6.2069.000

Park station for the measuring head of Professional VA/CVS instruments



1 PCS

6.2709.110

Stopper for the electrode opening on Professional VA /CVS instruments

Set with 1 stopper for electrode opening (e.g. for pH electrode or Pt1000 temperature sensor), 6 stoppers for openings (auxiliary solutions and/or automatic sample addition).



Optional accessories

Order no.	Description	
6.1454.130	Form seal made of FPM for 61414xxx for KF and VA instruments	
