



KF drying tube / SGJ 14 and straight

6.1403.010

Used for Karl Fischer water determination

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.1403.010

Qt.	Order no.	Description
1 PCS	6.2701.020	Cover
		Part of 6.1403.010 and 6.1609.000 drying tubes



Optional accessories

Order no.	Description	
6.1414.010	Titration vessel lid with 5 openings	
