



SPE measuring head insert for VA Computrace and 663 VA Stand

6.1412.020

Measuring head insert for operation with *screen-printed electrodes* (SPE) or the scTRACE Gold. With installation material.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 6.1412.020

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.1454.130

Form seal made of FPM for 61414xxx

for KF and VA instruments



1 PCS

6.1808.000

Connection bushing M6



1 PCS

6.2621.140

Hexagon key 2.5 mm



Optional accessories

Order no.	Description	
6.2709.130	Adapter B-14/15 - SPE measuring head insert (VA)	
	Adapter for use of the electrode shaft scTRACE Gold 6.1241.080 or the bismuth drop electrodes 6.0346.000 in the 884 Professional VA with SPE measuring head 6.1256.030 or in the 797 VA Computrace or 663 VA Stand with measuring head insert 6.1412.020.	

