



Titration vessel lid automation

6.1414.080

Titration vessel lid for external use with the 6.1415.250 titration vessel base in automation.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.1414.080

Qt.	Order no.	Description
1 PCS	6.1454.130	Form seal made of FPM for 61414xxx
		for KF and VA instruments



Optional accessories

Order no.	Description
6.2730.080	Screw nipple PCTFE. Used with 6.5613.000 ECO titration equipment with 722 Stirrer


