



KF titration vessel with lateral opening / 50-90 mL

6.1421.250

For Karl Fischer water determination Used together with spoon for solid substances 6.2702.000.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.1421.250

Qt.	Order no.	Description
1 PCS	6.2036.000	Mounting ring for titration vessel



Optional accessories

Order no.	Description	
6.2702.000	Weighing spoon with protective tube	
	For 6.1421.250 titration vessel. For KF water determination.	
