



Reaction vessel short for stability measurements

6.1429.040

Reaction vessel for stability measurements with Rancimat and PVC Thermomat. Set of 117 pieces.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.


Optional accessories

Order no.	Description	
6.1808.090	Thread adapter M8 outer / M6 inner Outer thread M8, inner thread M6, e.g. for mounting 6.2418.XXX Air tubes in 6.1429.XXX Reaction vessels of stability measuring instruments.	

6.2418.100

Air tube for oil measurements

Set of 117 pieces. In combination with 6.1429.040.



6.2418.120

Gas inlet tube for PVC measurements

Set of 117 pieces. In combination with 6.1429.040.

