



Stopper, PTFE, SGJ 29, for transfer set

6.1446.180

Stopper for titration and volumetric KF titration in the external titration vessel

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.1446.180

Qt.	Order no.	Description
2 PCS	6.1808.000	Connection bushing M6



1 PCS

6.2726.090

Antidiffusion valve

For antidiffusion buret tips 6.1541.XXX, 6.1543.130, 6.1543.140,
6.1543.150, 6.1543.160, 6.1543.200

