



Reagent Replacement Solvent Module stopper

6.1446.220

Stopper for reagent replacement in a coulometric titration cell with a Solvent Module.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.1446.220

Qt.	Order no.	Description
1 PCS	6.2713.000	PTFE sleeve SGJ14

For fat-free ground-joint connection with the standard ground-joint SGJ 14.

