



Robotic arm with Luer-Lock adapter for 786 Swing Head, right swinging

6.1462.090

Transfer robotic arm for liquid handling

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.1562.110	Pipetting tubing / 2 mL
	Pipetting tubing for liquid transfer



6.2816.100

Needle with vent

With luer lock, for transfer equipment 6.5619.000.

