



Robotic arm with holder for titration head and sensor for 786 Swing Head, right swinging

6.1462.170

Titration head holder with beaker sensor for 786 Swing Head

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.1458.010	Titration head, 6 x SGJ 14 and 3 x SGJ 9 openings. With 6xSGJ 14 and 3x SGJ 9 openings.



6.1458.020

Titration head / 4 x M10 connectors

With 4 M10 screw nipples, including stoppers.



6.1458.030

KF titration head / 4 x M10 connectors

With 4 M10 screw nipples, including stoppers.



6.1458.040

Titration Head, 3 x SGJ 14

Titration head for 120 mL beakers. With 3 openings SGJ 14 and 4 openings 6.4 mm, 2 dosing tips included. Suitable for rinsing and aspiration equipment.



6.1458.070

Titration head / 2 x SGJ 14

Titration head for 75 mL beakers. With 2 openings SGJ 14, 2 openings 6.4 mm, and 3 pre-mounted dosing tips. Suitable for rinsing and aspiration equipment.

