



Buret tip, tubing length 0.4 m

6.1541.010

With antidiffusion valve.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.1541.010

Qt.	Order no.	Description
1 PCS	6.1446.030	Ball stopper / B-14/(15) For the movable holding of buret tips in titration vessels



1 PCS

6.1543.200

Titration tip / M6 thread

Titration tip with antidiffusion valve.



1 PCS

6.1805.100

FEP tubing / M6 / 40 cm

With light and kink protection

