



Tip, length 280 mm

6.1543.130

With antidiffusion valve.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.1543.130

Qt.	Order no.	Description
1 PCS	6.2726.090	Antidiffusion valve
		For antidiffusion buret tips 6.1541.XXX, 6.1543.130, 6.1543.140, 6.1543.150, 6.1543.160, 6.1543.200

