



Titration tip / M6 thread / 440 mm

6.1543.180

Titration tip with antidiffusion valve.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.1543.180

Qt.	Order no.	Description
1 PCS	6.2726.090	Antidiffusion valve
		For antidiffusion buret tips 6.1541.XXX, 6.1543.130, 6.1543.140, 6.1543.150, 6.1543.160, 6.1543.200

