



3-way stopper with antidiffusion valve

6.1543.210

Stopper with three buret tips for connecting 3 x M6 tubing

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.1543.210

Qt.	Order no.	Description
3 PCS	6.2726.090	Antidiffusion valve For antidiffusion buret tips 6.1541.XXX, 6.1543.130, 6.1543.140, 6.1543.150, 6.1543.160, 6.1543.200

