



Pipetting equipment for Manual Dosing Controller

6.1562.060

Pipetting equipment for Manual Dosing Controller for 865 and 876 Dosimats.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.1562.140	Pipetting tubing 10 mL for 6.1562.060

A coiled white pipetting tube with a connector. The tube is white and flexible, with a black connector at one end. It is shown in a coiled position, with the connector pointing towards the bottom right.