




## Pipetting tubing / 10 mL

6.1562.100

Pipetting tubing for fully automatic transfer of sample solutions.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

### Optional accessories

| Order no.  | Description                                                                                                  |
|------------|--------------------------------------------------------------------------------------------------------------|
| 6.1823.010 | Guide shaft to pipetting tubing<br>For assembly in the transfer robotic arm on the Robotic Sample Processors |

