



Thermostat jacket for 806 Exchange Units

6.1563.040

For 6.3026.XXX Exchange Units

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.1563.040

Qt.	Order no.	Description
2 PCS	6.1808.050	Tubing adapter olive / M8 outer 1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of exchange units and stability measuring instruments.



