



Cylinder element ETFE 5 mL

6.1566.150

For dosing units used for dosing alkali reagents or hydrofluoric acid (HF) with Dosinos. Less precision and reproducibility is to be expected in comparison with glass cylinders.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.1566.150

Qt.	Order no.	Description
0 PCS	6.2803.000	Grease (1 oz.)
		Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX

