



Cylinder element glass 2 mL

6.1574.120

Glass cylinder element complete with cylinder, piston, cylinder base, valve disk, and centering tube. For 807 Dosing Unit 2 mL.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.1574.120

Qt.	Order no.	Description
0 PCS	6.2803.000	Grease (1 oz.) Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX

