



Cylinder element glass 50 mL

6.1574.250

Glass cylinder element complete with cylinder, piston, cylinder base, valve disk, and centering tube. For 807 Dosing Unit 50 mL.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.1574.250

Qt.	Order no.	Description
-----	-----------	-------------

0 PCS

6.2803.000

Grease (1 oz.)

Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX

