



## 807 Dosing Unit ETFE 50 mL

6.1575.250

807 Dosing Unit ETFE 50 mL for dosing alkali reagents or hydrofluoric acid (HF) with Dosinos. Less precision and reproducibility is to be expected in comparison with glass cylinders.

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

Scope of delivery 6.1575.250

| Qt. | Order no. | Description |
|-----|-----------|-------------|
|-----|-----------|-------------|

1 PCS

6.1446.040

Threaded stopper / M6



---

1 PCS

6.2803.010

Grease (2 g)

Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX



## Optional accessories

| Order no.  | Description                                                                                                                                                                                                          |                                                                                     |
|------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| 6.1574.250 | Cylinder element glass 50 mL<br><br>Glass cylinder element complete with cylinder, piston, cylinder base, valve disk, and centering tube. For 807 Dosing Unit 50 mL.                                                 |  |
| 6.1566.250 | Cylinder element ETFE 50 mL<br><br>For dosing units used for dosing alkali reagents or hydrofluoric acid (HF) with Dosinos. Less precision and reproducibility is to be expected in comparison with glass cylinders. |  |