

Cylinder element glass PCR 5 mL

6.1578.150

Glass PCR cylinder element complete with cylinder, piston, cylinder base, valve disk, and centering tube. For the 6.3033.150 dosing unit PCR 5 mL.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.1578.150

Qt.	Order no.	Description
0 PCS	6.2803.000	Grease (1 oz.) Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX

