



807 Dosing Unit 50 mL without accessories

6.1580.250

807 Dosing Unit with integrated data chip with 50 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.1580.250

Qt.	Order no.	Description
1 PCS	6.2803.010	Grease (2 g)
		Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX

