




Bottle attachment / S 40 / complete

6.1602.135

Enables direct dosing from reagent bottles with S 40 thread (Merck...).

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.1602.135

Qt.	Order no.	Description
1 PCS	6.1446.080	Threaded stopper / M8



1 PCS

6.1602.115

Bottle attachment / S 40 / green

For direct dosing from reagent bottles with S 40 thread (Merck...).



1 PCS

6.1805.080

FEP tubing / M6 / 25 cm

With light and kink protection.



1 PCS

6.1819.020

Tube

For M6 threaded holes. For 6.1602.1XX bottle attachments



1 PCS

6.2023.020

Clip for SGJ 14/15

Clip for SGJ 14/16

