



Bottle attachment GL 45 for eluents and reagents

6.1602.160

For eluent, reagent, and rinsing bottles, with connectors for the adsorber tube and the aspiration tubing.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.1602.160

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1602.105

Bottle attachment / GL 45 / green

For direct dosing from reagent bottles with GL 45 thread.



1 PCS

6.1446.080

Threaded stopper / M8

