



Bottle attachment GL 45 for level sensor

6.1626.000

This bottle attachment, fitted with the 6.2769.XXX sensors, results in the level sensors EMPTY and FULL, respectively, for bottles and canisters.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.1626.000

Qt.	Order no.	Description
2 PCS	6.1446.040	Threaded stopper / M6



1 PCS

6.1446.080

Threaded stopper / M8



1 PCS

6.1619.000

Adsorber tube for 807 Dosing Unit

