



Tubing connector for grounded buret tip

6.1808.030

Together with 6.1541.010 buret tip forms the complete grounded 6.1540.010 buret tip. For titrations in non-aqueous solutions.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 6.1808.030

Qt.	Order no.	Description
1 PCS	6.1808.010	Coupling piece 2 X M6
		2 M6 outer threads. Part of the screw connections for Exchange Units and 676 Sample Changer

