



Support stand with support rod

6.2001.010

For Dosimats and Titrinos

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.2001.010

Qt.	Order no.	Description
1 PCS	6.2016.050	Support rod / 300 mm



1 PCS

6.2621.100

3 mm hex key for IC Sample Processors



1 PCS

6.2621.070

5 mm hex key for IC Sample Processors

