



Stand plate for 912/913/914

6.2001.130

Stand plate for converting a mobile 912/913/914 pH/Conductometer into a laboratory meter.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.2001.130

Qt.	Order no.	Description
1 PCS	6.2013.010	Clamping ring

For support rods with a diameter of 10 mm.



1 PCS

6.2016.050

Support rod / 300 mm



1 PCS

6.2021.020

Electrode holder

Electrode holder for 4 electrodes and 2 buret tips



1 PCS

6.2621.070

5 mm hex key for IC Sample Processors

