



Sample rack, 48 x 11 mL, (Al)

6.2041.200

Sample rack for 889 IC Sample Center, suitable for 48 sample vessels. The following sample vessels can be used: 6.2743.100, 6.2743.107, 6.2743.110, 6.2743.117.

Material: Aluminum (Al)

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description
-----------	-------------

6.2743.100 PP Vials 0.7 mL to IC Sample Center, 1000 pieces
1000 sample vials.



6.2743.107 PP Vials 0.7 mL to IC Sample Center, 100 pieces
100 sample vials.



6.2743.110 PP Vials 0.3 mL to IC Sample Center, 1000 pieces
1000 sample vials.



6.2743.117 PP Vials 0.3 mL to IC Sample Center, 100 pieces
100 sample vials.



6.7303.100

Glass vial, 1.5 mL, 100 pieces


For use with the following sample changers: Liquid Autosampler CIC (TEI), MMS 5000 Autosampler (AJ) and 889 Sample Center.



6.7303.110

Septum for 67303100, 100 pieces

For use with the 67303100 Glass vial on the following sample changers: Liquid Autosampler CIC (TEI), MMS 5000 Autosampler (AJ), and 889 Sample Center.



6.2743.120

Septum to 6.2743.1XX Si/PTFE

1000 pieces



6.2743.127

Septum to 6.2743.1XX Si/PTFE

100 pieces

