




Sample rack, 28 x 50 mL, (PVC)

6.2041.490

Sample rack for 919 Autosampler plus, 858 Professional Sample Processor suitable for 28 sample vessels. The following sample vessels can be used: 6.2747.010.

Plastic: Polyvinyl chloride (PVC)

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
-----------	-------------

6.2747.010

Sample vessel 50 mL

With blue screw cap and mL graduation. Set of 25 pieces.



6.1458.080

Transfer head for 919 IC Autosampler plus for VA and CVS applications.

