



Sample rack, 36 x 6 mL, (Al)

6.2041.700

Sample rack for 774 Oven Sample Processor, suitable for 36 sample vessels. The following sample vessels can be used: 6.2419.000.

Material: Aluminum (Al)

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description
-----------	-------------

6.2419.000

Sample glass, 6 mL, 1000 pieces

Sample glass for up to 6 mL sample, suitable for 6.1448.050 septum,
1,000 pieces



6.1448.050

Aluminum septum caps / 1000 pieces

Including Al sealing.

