



Sample rack, 160 x 6 mL, (PVC)

6.2041.710

Sample rack for 815 Robotic Sample Processor, 855 Robotic Titrosampler "Basic" suitable for 160 sample vessels. The following sample vessels can be used: 6.2419.000.

Plastic: Polyvinyl chloride (PVC)

Below, the accessories are grouped into Scope of delivery and Optional accessories.


Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
-----------	-------------

6.1448.050 Aluminum septum caps / 1000 pieces
Including Al sealing.



6.1462.030 Robotic arm with transfer head for 786 Swing Head, left swinging
Transfer robotic arm for liquid handling purposes



6.1462.040 Robotic arm with transfer head for 786 Swing Head, right swinging
Transfer robotic arm for liquid handling purposes



6.1462.090 Robotic arm with Luer-Lock adapter for 786 Swing Head, right swinging
Transfer robotic arm for liquid handling



6.2419.000

Sample glass, 6 mL, 1000 pieces

Sample glass for up to 6 mL sample, suitable for 6.1448.050 septum,
1,000 pieces