



Sample rack, 28 x 250 mL, (PVC)

6.2041.820

Sample rack for 815 Robotic USB Sample Processor XL, 855 Robotic Titrosampler "Basic" suitable for 28 sample beakers. The following sample beakers can be used: 6.14332.320, 6.1432.323, 6.1453.220, 6.1453.250.

Plastic: Polyvinyl chloride (PVC)

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
-----------	-------------

6.1432.320 Sample beaker, clear glass, 250 mL, 1 piece

Sample beaker, clear glass, 250 mL, suitable for sample racks:
6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820,
6.2068.020 and 6.2068.030, 1 piece



6.1453.220 Sample beaker, plastic (PP), 200 mL, 1 piece

Sample beaker, plastic (polypropylene, PP), 200 mL, suitable for sample
racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820,
6.2068.020, and 6.2068.030, 1 piece



6.1453.250 Sample beaker, plastic (PP), 250 mL, 1 piece

Sample beaker, plastic (polypropylene, PP), 250 mL, suitable for sample
racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820,
6.2068.020, and 6.2068.030, 1 piece



6.1458.010 Titration head, 6 x SGJ 14 and 3 x SGJ 9 openings.

With 6xSGJ 14 and 3x SGJ 9 openings.



6.1458.040

Titration Head, 3 x SGJ 14

Titration head for 120 mL beakers. With 3 openings SGJ 14 and 4 openings 6.4 mm, 2 dosing tips included. Suitable for rinsing and aspiration equipment.



6.1462.060

Robotic arm with holder for titration head, left swinging

Titration head holder to 786 Swing Head



6.1462.070

Robotic arm with holder for titration head, right swinging

Titration head holder to 786 Swing Head



6.1462.160

Robotic arm with holder for titration head and sensor for 786 Swing Head, left swinging

Titration head holder with beaker sensor for 786 Swing Head



6.1462.170

Robotic arm with holder for titration head and sensor for 786 Swing
Head, right swinging

Titration head holder with beaker sensor for 786 Swing Head

