




Bottle holder to Extension Module - Liquid Handling

6.2057.200

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.2057.200

Qt.	Order no.	Description
1 PCS	6.2621.030	Hexagon key 4 mm



Optional accessories

Order no.	Description	
6.1562.160	Transfer tubing 15 mL, 2 x M6 For liquid handling in ion chromatography. To avoid contamination of the Dosino.	
6.1608.100	PE bottle 50 mL Bottle for Extension Module Liquid Handling.	
6.2762.000	Mixing vessel 50 mL to Extension Module - Liquid Handling	