



Fastening plate for Dosino holder

6.2058.120

Fastening plate for flexible connection of the Dosino holder at the 862 Compact Titrosampler. Thanks to the arbitrary mounting height it is easy to connect reagent bottles of different suppliers.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 6.2058.120

Qt.	Order no.	Description
1 PCS	6.2621.130	Hexagon key 2 mm
	2 mm.	

