



Sample rack, 24 x 120 mL, (PVC)

6.2068.020

Sample rack for 815 Robotic Soliprep, suitable for 24 sample beakers. The following sample beakers can be used: 6.1459.300, 6.1459.307.

Plastic: Polyvinyl chloride (PVC)

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
-----------	-------------

6.1459.300 Sample beaker, plastic (PP), 120 mL, 250 pieces

Sample beaker, plastic (polypropylene, PP), 120 mL, suitable for sample racks: 6.2041.470, 6.2041.840, 6.2041.850, 6.2068.010, 6.2068.020, 6.2068.030, and 6.9920.141, 250 pieces




6.2067.000 Insert for 30x11 mL for Soliprep sample rack

Sample rack insert for 30x11 mL (6.2743.050)



6.2067.010 Insert for 24xSyringe needles for Soliprep sample rack

Sample rack insert for 24xsyringe needles (luer connection)



6.2067.020 Insert for 12x30 mm syringe filters for Soliprep Sample rack

Sample rack insert for 12x30 mm syringe filters (luer connection)



6.1459.307

Sample beaker, plastic (PP), 120 mL, 1 piece

Sample beaker, plastic (polypropylene, PP), 120 mL, suitable for sample racks: 6.2041.470, 6.2041.840, 6.2041.850, 6.2068.010, 6.2068.020, 6.2068.030, and 6.9920.141, 1 piece

