



Sample rack, 10 x 250 mL, (PVC)

6.2068.030

Sample rack for 815 Robotic Soliprep, suitable for 10 sample beakers. The following sample beakers can be used: 6.1432.320, 6.1432.323, 6.1453.220, 6.1453.250.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.1432.320	Sample beaker, clear glass, 250 mL, 1 piece
	Sample beaker, clear glass, 250 mL, suitable for sample racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820, 6.2068.020 and 6.2068.030, 1 piece



6.2067.000 Insert for 30x11 mL for Soliprep sample rack
Sample rack insert for 30x11 mL (6.2743.050)




6.2067.020 Insert for 12x30 mm syringe filters for Soliprep Sample rack
Sample rack insert for 12x30 mm syringe filters (luer connection)



6.2067.030 Insert 30x11.6 mm for Soliprep Sample rack
Sample rack insert for 30 x 11.6 mm vials. The upper part for fixing the vials is already mounted on the Soliprep sample rack (6.2068.030).



6.2067.040 Insert 30x syringe needles for Soliprep sample rack
Sample rack insert for 30 x syringe needles (luer connection)



6.1432.323 Sample beaker, amber glass, 250 mL, 1 piece

Sample beaker, amber glass, 250 mL, suitable for sample racks:
6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820,
6.2068.020, and 6.2068.030, 1 piece



6.1453.220 Sample beaker, plastic (PP), 200 mL, 1 piece

Sample beaker, plastic (polypropylene, PP), 200 mL, suitable for sample
racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820,
6.2068.020, and 6.2068.030, 1 piece



6.1453.250 Sample beaker, plastic (PP), 250 mL, 1 piece

Sample beaker, plastic (polypropylene, PP), 250 mL, suitable for sample
racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820,
6.2068.020, and 6.2068.030, 1 piece

