



RS-232 cable 1 m, 25-pin D – sub 25

6.2125.090

RS-232 extension cable (1:1). Specially shielded.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.2125.120	Adapter sub D 25-pin/5 sockets Distributes the «Remote lines to 5 x 25-pin sockets. Used with pH Meters, Titrinos, KF Coulometer, Sample Changers, IC systems, etc.

