



Buffer solution pH 4 (500 mL)

6.2307.100

Ready-to-use buffer solution pH = 4 (25 °C), colored, with closure seal

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.2307.230	Buffer solutions pH 4, 7 and 9 Mixed buffer solutions pH 4.00/7.00/9.00 (25 °C) in single use sachets, colourless, box of 3 x 10 x 30 mL



6.2307.210 Buffer solution pH 7.00 (30 x 30 mL)
Buffer solution pH 7.00 (25 °C) in single use sachets, colorless



6.2307.220 Buffer solution pH 9.00 (30 x 30 mL)
Buffer solution pH 9.00 (25 °C) in single use sachets, colorless



6.2307.110 Buffer solution pH 7 (500 mL)
Ready-to-use buffer solution pH = 7 (25 °C), colored, with closure seal



6.2307.120 Buffer solution pH 9 (500 mL)
Ready-to-use buffer solution pH = 9 (25 °C), colorless, with closure seal



6.2325.000

pHit kit

Maintenance kit for pH electrodes

The kit contains:

- 50 mL cleaning solution
- 50 mL 3M KCl solution
- 50 mL storage solution
- 2 Storage vessels
- Instructions for use



6.2307.200

Buffer solution pH 4.00 (30 x 30 mL)

Buffer solution pH 4.00 (25 °C) in single use sachets, colorless

