



## Capillary cutter

6.2621.080

Cutting tool for simple cutting of plastic capillaries, 1/16" and 1/8", including 5 replacement blades.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

### Optional accessories

| Order no.  | Description                                                              |
|------------|--------------------------------------------------------------------------|
| 6.1803.040 | PTFE capillary 0.5 mm i.d. / 1 m<br>Capillary for sample handling in IC. |



6.1803.080

PTFE capillary 0.97 mm i.d. / 1 m

Capillary for filtration, for IC Filtration Sample Processor, IC Liquid Handling Filtration Unit



---

6.1831.010

PEEK capillary 0.25 mm i.d. / 3 m

For all IC components.



---

6.1831.030

PEEK capillary 0.75 mm i.d. / 3 m

For all IC components.



---

6.1831.060

PEEK capillary 0.5 mm i.d., 1 m

For inline dialysis and for Sample Changers. For IC Dialysis Unit, IC Sample Processor, Compact IC Autosampler, IC Filtration Sample Processor, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit

