



Drip pan for 789/815/855

6.2711.070

Drip pan for 789, 815, 855 Robotic Sample Processors XL

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.2711.070

Qt.	Order no.	Description
1 PCS	6.1801.130	PVC tubing / 8 mm / 12 mm / 2 m



Optional accessories

Order no.	Description	
6.2058.030	Head extension 10 mm for Sample Processors Head extension for use with the waterbath rack.	
