



Filtration membrane reg. cellulose, 0.45 µm, 50 pieces

6.2714.050

Nominal pore size: 0.45 µm. Set of 50 pieces. For inline filtration.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories


Order no.	Description
6.02729.110	Ultrafiltration cell 2

For the Inline ultrafiltration of samples together with the 858 Professional Sample Processor / 919 IC Autosampler plus. As an alternative to the time-consuming manual filtration with a syringe filter.



6.05330.010 IC equipment: Inline Ultrafiltration 2 - pull mode


Accessory set for Inline Ultrafiltration 2 - pull mode. For use with the 858 Professional Sample Processor/ 919 IC Autosampler plus.



6.05330.110 IC equipment: Inline Ultrafiltration 2 - push mode

Accessory set for Inline Ultrafiltration 2 - push mode. For use with the 858 Professional Sample Processor/ 919 IC Autosampler plus.

Recommended set for the combination with Inline Dilution (MIDT).



6.05330.210 IC equipment: Inline Ultrafiltration 2 - MiPT

Accessory set to combine Inline Ultrafiltration 2 and MiPT. For use with the 858 Professional Sample Processor.



6.2831.010 Plastic tweezers

For dialysis and filtration membranes

