




Dialysis cell

6.2729.100

With 4 x 6.2744.060 Blind stoppers. For inline dialysis

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.2729.100

Qt.	Order no.	Description
4 PCS	6.2744.060	Threaded stopper
		For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Optional accessories

Order no.	Description
6.2744.000	Pressure screw PVDF With UNF 10/32 connection. For IC instruments, inline dialysis (connection of PTFE capillaries)

