




Safety shield for Robotic Soliprep

6.2751.150

Safety shield for the Robotic Soliprep systems Filtration and Flexible, which pick up and then discard various tools at Tower 1, e.g. syringe needles. Reaching into the movement radius of the robotic arm is prevented and at the same time protection is offered against inadvertently escaping solutions.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.2751.150

Qt.	Order no.	Description
-----	-----------	-------------

