



FULL sensor for level sensor

6.2769.000

This sensor, mounted in the 6.1626.000 bottle attachment, results in the level sensor FULL for 1 L and 2 L bottles and for canisters. Included are two sensor rods, one inlet tubing and one holder for fixing the tubing to one of the sensor rods.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.2769.000

Qt.	Order no.	Description
1 PCS	6.1829.010	FEP aspiration tube M6, 0.25 m With M 6 thread. For Dosing Unit

