



Polishing set for Ultra Trace Electrode Tip

6.2802.020

Polishing set to renew the surface of the Ultra Trace Electrode tip 6.1204.180

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.2802.020

Qt.	Order no.	Description
1 PCS	6.2802.030	Polishing tape to Polishing Set for Ultra Trace Electrode
		Polishing tape to Polishing Set for Ultra Trace Electrode. 12 µm particle size. used in the polishing cylinder. Set of 50

