



# Syringe 10 mL with Luer connection

6.2816.020

For various applications in IC and VA

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

## Optional accessories

| Order no.  | Description                                                                     |
|------------|---------------------------------------------------------------------------------|
| 6.2744.120 | Luer F pressure screw adapter<br>Connection between pressure screw and syringe. |



6.2816.030

Needle with Luer connector

Needle for filling in mercury into the Multi-Mode Electrode and for sample addition in Karl Fischer titration.



---

6.2816.040

Purge tubing

With PTFE tubing and Luer connection. For syringes. For aspirating eluents.

