



Standard pump head PEEK

6.2824.110

Standard pump head for intelligent IC instruments

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.2824.110

Qt.	Order no.	Description
2 PCS	6.2741.020	Piston seal (orange) For all standard pump heads



1 PCS

6.2824.170

Inlet valve PEEK

Metal-free check valve for the inlet side of the IC pump heads of the devices 940, 930, 925, 883, 882, 881, and 850.



1 PCS

6.2824.160

Outlet valve PEEK

Metal-free check valve for the outlet side of the IC pump heads of the devices 940, 930, 925, 883, 882, 881, and 850.

