



Water bath for Robotic Sample Processors and Robotic Titrosampler

6.2840.000

Water bath rack for Robotic Sample Processors and Titrosampler.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.2840.000

Qt.	Order no.	Description
2 PCS	6.1801.130	PVC tubing / 8 mm / 12 mm / 2 m



Optional accessories

Order no.	Description	
6.2041.900	Sample rack insert, 54 x 75 mL, for water bath, (Al) Sample rack insert for 815 Robotic USB Sample Processor XL, 855 Robotic Titrosampler "Basic", suitable for 54 sample beakers. The following sample beakers can be used: 6.1432.210, 6.1459.400. Material: Aluminum (Al)	
6.2058.030	Head extension 10 mm for Sample Processors Head extension for use with the waterbath rack.	